

Lipid Nanoparticles Deliver Nucleic Acids, Including CRISPR Components, in Primary Neurons

Precision NanoSystems Inc., Vancouver, BC, Canada, Genome Ctr., UC Davis, Davis, CA, Department of Neurology and Stem Cell Program, UC Davis Med. Ctr., Sacramento, CA


Grace Tharmarajah¹, Anitha Thomas¹, Rebecca De Souza¹, Ian Backstrom¹, Eric Ouellet¹, Shyam Garg¹, Keara Marshall¹, Shannon Chang¹, Timothy Leaver¹, Andre Wild¹, Peter Deng², Kyle D. Fink³, R. James Taylor¹ and Euan Ramsay¹

Introduction

Advances in the gene editing arena, specifically with CRISPR-Cas9, has pushed the demand for efficiently delivering payloads even further. Of the tools available, developments in the field of lipid nanoparticles (LNPs) has allowed for the reliable and efficient delivery of CRISPR components, both in research and clinical settings. Here, we bridge that gap by describing the development of an LNP delivery system for CRISPR components, robustly manufactured with clinical-grade materials using microfluidic technology at scales for in vitro applications, with the potential to move into animal models. We describe the use of lipid-based nanoparticles for highly efficient encapsulation and delivery of payloads, such as siRNA, mRNA, and plasmid. We show that representative small RNAs, mRNAs and plasmids can successfully be delivered to primary neurons. LNPs manufactured to encapsulate various nucleic acids can do so with high efficiency, encapsulating more than 95% of the payload, minimizing payload loss. Transfection efficiency of the LNPs is variable based on payload. The biological endpoint assays used to determine the accessibility of the payloads delivered varies for siRNA, mRNA, and plasmid. Using doses of 1 µg per mL of media, we achieved >90% knockdown with siRNA delivery, >90% of the primary neurons are GFP+ with GFP mRNA delivery and 60% of the primary neurons are GFP+ with GFP plasmid delivery. The LNPs are well tolerated, such that 5X the required dose have no observable cytotoxicity. We demonstrated delivery of gRNA to Cas9 expressing cells to achieve up to 80% editing. We demonstrated co-transfection of Cas9 mRNA and gRNAs in separate nanoparticles to achieve editing in wild-type cells. We show that the LNPs can also be used to deliver payloads into various regions of the animal brain. These validation studies provide suitable insights in establishing strategies for efficiently delivering CRISPR components into primary cultures.

Gene Delivery via Lipid Nanoparticles


Lipid Nanoparticles (LNPs) are a platform that can be used to deliver nucleic acids to cells. LNPs mimic low density lipoproteins (LDLs), which are taken up by an endogenous pathway. LNPs are pH sensitive, and designed to release their payload into the cytoplasm.




The Microfluidic Platform

The NanoAssemblr™ platform uses microfluidics to precisely control the manufacturing of nanoparticles from µL to mL scales.


in vitro scale




in vivo scale



Proof of Concept: mRNA Delivery in Rat Primary Neurons




MAP2 positive neurons (red) show expression of GFP (green) after treatment with Neuro9 GFP mRNA-LNP in the presence of 5 µg/mL ApoE. Cells were stained with anti-MAP2 primary antibody/alexa fluor-594 conjugated secondary and anti-GFP primary antibody/alexa fluor-488 conjugated secondary antibody. Nuclei are stained with DAPI (blue)



The percentage of GFP+ neurons was dose dependent, with 1 µg/mL having the highest percentage.

Flow cytometry analysis after 48 h treatment with Neuro9 GFP mRNA-LNP in the presence of 5 µg/mL ApoE in rat primary neurons, DIV 7, show >95% of neurons had taken up the LNP (calculated using the fluorescent probe DiD incorporated within the nanoparticle) even when treated at different doses of GFP mRNA LNP (Data not shown).



Cell viability was assayed by the PrestoBlue® Cell Viability Assay. No toxicity was observed at any dose employed, n=6 for each treatment One-way ANOVA with Dunnett's multiple comparison test.

The size and polydispersity index (PDI) of the nanoparticles encapsulating different mRNAs are comparable. GFP mRNA is 996 nucleotides long, while Cas9 mRNA is 4479 nucleotides in length. The particle size and PDI is not significantly different despite using mRNAs of different lengths, student's t-test.

Methods

Primary Rat Cortical Neuron Culture

E18 rat cortical tissue was purchased from Brainbits, LLC. The cortices were removed from the imaging media and transported with 0.25% trypsin-EDTA (ThermoFisher). The trypsin-EDTA was then inactivated by washing the tissue with 10% FBS in DMEM (ThermoFisher), followed by DNEM alone. The tissue was triturated in Neurobasal Media (NeuroCult™ Neuronal Base Medium (StemCell) with NeuroCult™-SM1 Neuronal supplement (StemCell) and L-Glutamine (StemCell) supplemented with L-Glutamic Acid (Sigma), and then passed through a 40µm cell strainer to form a single cell suspension. This suspension was then counted and plated on PDLC coated culture plates or glass cover slips at a density of 4x10⁴ cells/cm². Half of the media was changed with neuronal media without L-Glutamic Acid every 3-4 days.

Treatment of Primary Rat Cortical Neurons with LNPs

After 7 days in culture (DIV7) the 5µg/mL of ApoE was added following the specified LNP at the specified dose. The neurons were then incubated for 48h after which they were harvested for endpoint assessment.

Flow Cytometry

Following treatment and incubation with the specified LNP as noted above the neurons were washed with PBS and then trypsinized with 0.25% trypsin-EDTA to detach them from the culture plate. Trypsinization was inactivated by 3% FBS in PBS and the neurons were triturated to form a single cell suspension. The neurons were then pelleted and washed with PBS before being re-suspended in suspension buffer (BD Biosciences) and Propidium Iodide (BD Biosciences) was added to stain for apoptotic cells. The suspended, stained neurons were then passed through a 40µm cell strainer and then assayed using a BD Celesta Flow Cytometer.

Viability Assay

The standard protocol included with the PrestoBlue® Cell Viability Reagent (ThermoFisher), following the manufacturer's instructions. Values were corrected to wells containing media without neurons.

Immunocytochemistry

For this assessment the neurons were plated on PDLC coated coverslips, following cultivation, treatment and incubation, the neurons were washed with PBS and fixed with 4% PFA. The neurons were permeabilized using 0.1% Triton-X and then blocked with NGS followed by incubation with a MAP2 antibody (Sigma M4403) for 1 hour. The neurons were then blocked in NGS followed by incubation with a GFP antibody (AbCam ab19970) or Cas9 antibody (Diagene C15200229) in NGS overnight. Secondary antibodies were then added, AlexaFluor-594 for MAP2 and Cas9 and AlexaFluor-488 for GFP. The coverslips were mounted on glass slides using ProLong® Diamond Antifade Mountant (ThermoFisher) and imaged on a confocal microscope. Confirmation of the selective binding of the primary and secondary antibodies was provided through no primary controls (data not shown).


In Vivo Plasmid Transfection

VAC128495 mice at 8 months of age were intercranially injected with LNP-encapsulating plasmid encoding a protein with an eGFP reporter. Mice were anesthetized with 1% Isoflurane. The heads were prepared and a single craniotomy was made ~0.5mm anterior and ~1.75mm lateral to bregma. A Hamilton syringe was loaded with 5µL of plasmid-LNP and lowered ~2.5mm from the dura. 5µL of plasmid-LNP was injected at rate of 0.5 µL/minute. Needle was left in place for 5 min before being retracted over course of 2 minutes. Skull was sealed with bone wax and sutured. Mice were observed for recovery until ambulatory. 48 hrs following injection mice were intercranially perfused with PBS and formalin. Brains were extracted and cryoprotected in 30% sucrose before flash freezing in isopentane. Brains were sectioned using a cryostat at 30 µm thick serial sections every 4th section being collected for analysis. Brains were mounted and cross-labeled with Hoechst (1:1000) for 15 min. Slices were imaged (exposure 82-9000) at 4x and stitched together using Keyence Viewing Software.


TBEI Cleavage Assay

Following treatment and incubation of the Cas9-HEK293 cells or primary rat neurons, genomic DNA extraction and the subsequent TBEI cleavage assay was performed using the Alt-R® Genome Editing Detection Kit (IDT), following the manufacturer's instructions.

Proof of Concept: Plasmid Delivery in Rat Primary Neurons




MAP2 positive neurons (red) show expression of GFP (green) after treatment with Neuro9 GFP plasmid-LNP in the presence of 5 µg/mL ApoE. Cells were stained with anti-MAP2 primary antibody/alexa fluor-594 conjugated secondary and anti-GFP primary antibody/alexa fluor-488 conjugated secondary antibody. Nuclei are stained with DAPI (blue)



The percentage of GFP+ neurons was dose dependent, with 0.6µg/mL having the highest percentage.


Flow cytometry analysis after 48 h treatment with Neuro9 GFP plasmid-LNP in the presence of 5 µg/mL ApoE in rat primary neurons, DIV 7, show >95% of neurons showed uptake of LNPs (calculated using the fluorescent probe DiD incorporated within the nanoparticle) even when treated at different doses of GFP plasmid (data not shown).



Cell viability was assayed by the PrestoBlue® Cell Viability Assay. No toxicity was observed at any dose employed, n=6 for each treatment, one-way ANOVA with Dunnett's multiple comparison test.


The size and PDI of the nanoparticles encapsulating different plasmids are comparable. The particle size and PDI is not significantly different despite using plasmids of different size, student's t-test.

Plasmid Delivery in Vivo




GFP expression (green) 48H following striatal injection of GFP plasmid-LNP in 10 month old, FVB mice. Injection of 5µL of a 3mg/mL GFP plasmid-LNP Formulation. Co-localization of nuclear staining with DAPI (blue), LNP uptake using DiD (red), and GFP expression (green) in the same image. Both LNP uptake and GFP expression co-localize both within the area of injection, the striatum, and have spread to additional surrounding regions including parts of the ventricular system and the corpus callosum.

Delivery of Cas9 mRNA in Rat Primary Neurons



Primary rat neurons show Cas9 expression (Red) after treatment with Cas9 mRNA LNPs in the presence of 5µg/mL of ApoE. 48H after treatment the cells were stained for DAPI (blue) to identify nuclei, as well as Cas9 primary antibody/alexa fluor-594 conjugated secondary. At this time point the majority of the Cas9 protein has translocated to the nucleus indicating that delivery of guide RNA for the Cas9 protein must occur before this time point to ensure complexation and co-translocation to the nucleus.

Cas9 Mediated Editing in Cas9-HEK293 & Rat Primary Neurons



Cas9 expressing HEK293 cells were treated with three different formulations of hprt guide RNA LNPs. 48H following treatment the cells were harvested and efficiency of editing was determined using the TBEI cleavage assay.

Rat primary neurons were treated with 1µg/mL of Cas9 mRNA LNPs and 3µg/mL of guide LNPs in the presence of 1µg/mL of ApoE. 48H after treatment cells were harvested and efficiency of editing was determined using the TBEI cleavage assay.

Addition of guide RNA to Cas9 expressing cells produces robust editing in HEK293 cells. In wild-type neuronal cells further titration is necessary to ensure the correct timing for optimal editing

Conclusion

The NanoAssemblr™ platform is capable of creating highly reproducible LNPs carrying siRNA, mRNA or plasmid payloads.

These LNPs have high transfection efficiencies (>95%) in primary, rat, cortical neurons and allow for efficient delivery of their associated payloads.

The LNPs were not toxic to these primary neurons even at high dosages (50x)

LNP-mediated CRISPR-Cas9 Delivery

LNPs are versatile, and can be used to deliver Cas9 via mRNA, plasmid, or a pre-formed ribonucleoprotein complex. One or more guide strands can be delivered as well.

Our results validate that these lipid nanoparticles can be used as an effective delivery system of CRISPR-Cas9 components to mediated gene editing.

