# Low Volume Production Of Nanoparticles That Are Effective Transfection Systems In iPSC-derived Cells, Immune Cells and Other Primary Cell Cultures

Andrew Brown, Chelsea Cayabyab, Rebecca A. G. De Souza, Tara Fernandez, Grace Tharmarajah, Anitha Thomas, Samuel Clarke and Euan Ramsay


# Purpose and Objectives

- Nucleic acid therapies permit access to previously undruggable pathways to treat or potentially cure a wide range of diseases, and lipid nanoparticles (LNPs) are the most clinically advanced non-viral delivery systems
- Significant need exists for robust manufacturing of LNPs at all scales, particularly at very low volumes to enable efficient screening of lipid components and nucleic acid payloads at earlier stages of drug development
- Here, we demonstrate the use of a low volume microfluidic device that provides consistent, controlled formulation conditions for encapsulating just tens of micrograms of mRNA and plasmids into LNPs
- LNP materials were systematically screened for size, polydispersity index, and activity in clinically relevant induced pluripotent stem cell (iPSC)-derived neuroprogenitor cells (NPCs) and cortical neurons

# Methods

# Microfludic Nanoparticle Production


The NanoAssemblr<sup>™</sup> platform uses microfluidics to precisely control the manufacturing of nanoparticles from microlitre to litre scales.



- mediate endosomal disruption and release of nucleic acid
- LNPs were tested for expression in clinically relevant iPSC-derived NPCs and cortical neurons

# Screening microgram mRNA formulations

# Highest PEG concentration reduces encapsulation efficiency



### Formulation Details


Lipid Composition	PNI-ILa:DSPC:Cholesterol:PEG-DMG 50:10:40-x:x   $x = 1.5$ , 3, 5 mol%	
Init. Lipid Conc.	50 mM	
mRNA	GFP (996 nt), 0.025 mg	
N/P Ratio	4	
Organic Solvent	Ethanol	
Aqueous Phase	0.695 mg/mL mRNA in 100mM Sodium Acetate	
Instrument Setting	3x	

Downstream Processing 4x Dilution in nuclease-free PBS

### Further Details

Size measured by DLS and Encapsulation Efficiency by Ribogreen assay. Values represent the mean of triplicate formulations. Error bars are the standard deviation. Means grouped by a horizontal line are not significantly different (Tukey's multiple comparison test, P>0.05).

# Ionizable cationic lipid produces LNPs with lower size and PDI



Ionizable/cationic lipid


# Formulation Details

Lipid Composition	CL:DSPC:Cholesterol:PEG-DMG 50:10:48.5:1.5 mol% CL = cationic lipid as labeled on x-axis
Init. Lipid Conc.	50 mM
mRNA	GFP (996 nt), 0.025 mg
N/P Ratio	4
Organic Solvent	Ethanol
Aqueous Phase	0.695 mg/mL mRNA in 100mM Sodium Acetate
Instrument Setting	3x
Downstream Processing	4x Dilution in nuclease-free PBS

# Further Details

Size measured by DLS and Encapsulation Efficiency by Ribogreen assay. Values represent the mean of triplicate formulations. Error bars are the standard deviation and are not drawn when error is less than the size of the symbol. Means grouped by a horizontal line are not significantly different (Tukey's multiple comparison test, P>0.05).





Lipid Composition	PNI-ILa:DSI 50:10:38.5:	PC:Cholester :1.5 mol%	ol:PEG-DMG
Init. Lipid Conc.		50 mM	
N/P Ratio	4	6	8
mRNA		GFP (996nt)	
Mass (mg)	0.025	0.0167	0.0125
Conc. (mg/mL)	0.695	0.464	0.347
Aqueous Phase	mRNA in 10	mRNA in 100 mM Sodium Acetate	
Organic Solvent	Ethanol	Ethanol	
Instrument Setting	3x		

Downstream Processing 4x Dilution in nuclease-free PBS

# Further Details

Size measured by DLS and Encapsulation Efficiency by Ribogreen assay. Values represent the mean of triplicate formulations. Error bars are the standard deviation and are not drawn when error is less than the size of the symbol. Means grouped by a horizontal line are not significantly different (Tukey's multiple comparison test, P>0.05).




# Case Study: Optimizing LNP formulations for plasmid expression in iPSCs

# Formulations with similar properties perform differently in iPSC-derrived NPCs undergoing differentiation: Biological activity is a necessary screen


- Physical properties similar between formulations with different lipids and N/P ratios
- A panel of formulations containing PNI-ILa with different helper lipids was created at two N/P ratios
- Formulations 2, 3, and 4 have similar properties
- Panel was created rapidly with the Spark at µg quantities with encapsulation efficiency > 80%

### Further Details

Size measured by DLS and Encapsulation Efficiency by Ribogreen assay. Values represent the mean of duplicate formulations. Error bars are the standard deviation and are not drawn where error is smaller than the symbol.








# Systematic screening of plasmid LNP formulations in mature human iPSC-derived cortical neurons: N/P ratio had a profound impact

Hours Post Plasmid LNP Treatment



# Conclusions

- Microliter formulations containing microgram quantities of mRNA were rapidly produced using the NanoAssemblr Spark microfluidic system to systematically screen compositions against properties and activity
- PEG-lipid content  $\geq$  4 mol% reduced encapsulation efficiency; PNI-ILa produced LNPs with more favourable size and PDI than other cationic lipids; N/P ratios > 6 did not appreciably affect properties
- A panel of 4 formulations formulated at N/P 6 and N/P 4 were rapidly developed for testing in iPSC-derived neurons



This poster is interactve! Download the free "layar" app on your smartphone or tablet and use it to scan this poster. Enriched features are indicated with this logo.


Precision NanoSystems Inc, Vancouver, BC, Canada

# **B** Formulations 3 & 4 Exhibit the Least Impact on Neurite Health and

#### presentative Image Phase - Formulation No.4 @ 78h

- Formulation composition affects neuronal development: formulations 1 and 2 exhibit large decrease in neurite length that is not recovered over time
- Reporter gene expression peaks ~78h after treatment and appreciable expression is still present 14d following treatment
- The same formulations that do not impact neurite length also display high GFP expression, with approximately 80% of cell bodies expressing GFP by 78 hours

# **Treatment Details**



Time lapse video of iPSC-NPCs growing and expressing GFP over 21 day period: Scan with Layar app.

Cell Type	iPSC-Derived NPCs
N/P Ratio	6
Payload	GFP Plasmid
Treatment Point	Day 1 of differentiation protocol
Dose	0.6 µg DNA /mL media
Аро-Е	1 µg/mL media
<b>Biological Assay</b>	Live-cell imaging: Phase: Neurite Length Fluorescence: GFP Expr.
Capture Frequency	6 Hours
Total Observed Period	21 Days

#### • The N/P ratio is important to formulation efficacy for plasmid delivery to human iPSC-derived cortical neurons.

• Different N/P ratios were tested to identify candidate formulations with optimal performance in this therapeutically

# reatment Details

Cell Type	iPSC-Derived NPCs
N/P Ratio	4
Payload	GFP Plasmid
Treatment Point	Day 14 of differentiation
Dose	0.6 µg DNA /mL media
Apo-E	1 µg/mL media
Biological Assay	Live-cell imaging: Phase: Neurite Length Fluorescence: GFP Expr.
Capture Frequency	6 Hours
Total Observed Period	7 Days

- On iPSC-derived NPCs, formulations 3 and 4 outperformed 1 and 2 in terms of cell tolerance and reporter gene expression
- On iPSC-derived cortical neurons, all formulations at N/P 4 were better tolerated than N/P 6, and formulations 3 and 4 at N/P = 4 were most effective in inducing reporter gene expression
- Despite similar physical properties, formulations had very different in vitro efficacy
- The Spark is ideal for efficient screening of formulations with valuable payloads in clinically relevant cell types

