

Novel Lipid Nanoparticle Delivery Reagent and Rapid Manufacturing Workflow to Accelerate Preclinical Development of RNA Vaccines

Edward Grahame^{1,2}, Natalie Orr^{1,2}, Suraj Abraham¹, Sitalakshmi Thampatty¹, Rita Zhao¹, Angela Zhang¹, Andy Geall¹, Anitha Thomas¹, Samuel Clarke¹ and Yvonne Perrie²

¹Precision Nanosystems Inc., Vancouver, Canada


²University of Strathclyde, Glasgow, Scotland

Contact us at:
Precision NanoSystems Inc, Vancouver, BC, Canada
info@precision-nano.com


Purpose

- The biggest vaccine campaign in history is ongoing against COVID-19, with RNA vaccines playing a key role.
- Due to the rapid degradation of RNA and low cellular uptake, an essential element of RNA vaccines are lipid nanoparticle (LNP) delivery systems.
- Limited access to ionizable lipids and LNP compositions, and the difficulty in scale-up production of RNA-LNPs remain challenges in the field.
- In this study, we aim to highlight that commercially available LNP reagent mix, GenVoy-ILM, is an accessible and easy-to-use LNP formulation that allows for rapid preclinical development of RNA vaccines.

Objectives




GenVoy-ILM LNPs can encapsulate both mRNA and saRNA while retaining the CQAs required for an RNA vaccine.





GenVoy-ILM is an Effective *In Vitro* Delivery Vehicle for mRNA and saRNA

mRNA Expression (A): Dose-response of GenVoy mRNA treated HEK293 cells showing live/dead cell proportions, transfection efficiency (GFP%(live)) and representative histogram. **saRNA expression (B):** Dose-response of GenVoy saRNA treated HEK293 cells showing live/dead cell proportions, transfection efficiency (GFP%(live)) and representative histogram. N=3 ± SD




GenVoy-ILM LNPs induce an immune response against the target antigen


LNP vaccines were formulated with saRNA encoding for the full-length SARS-CoV-2 spike protein. (A) ELISA analysis of total anti-spike serum IgG levels demonstrated that GenVoy-ILM vaccine produced antibodies against the target antigen. (B) ELISA analysis of the anti-spike IgG2a and IgG1 isotypes demonstrated that the IgG response to the GenVoy-ILM vaccine was skewed towards the IgG2a isotype which is indicative of a Th1 immune response



Genetic Medicine Toolkit




Methods




RNA-LNP characterization and *in vitro* activity: RNA-LNP size, PDI and zeta potential were determined using DLS and ELS (Malvern Zetasizer Ultra). The encapsulation efficiency (EE%) of the RNA was determined using Ribogreen™ reagent. In vitro expression/potency and viability were determined by flow cytometry (Attune NxT). HEK293 cells were treated with eGFP mRNA and saRNA LNPs in the presence of ApoE (1µg/mL) and fluorescence was determined.

***In vivo* expression and immunogenicity:** Female BALB/c mice (n=5) were injected IM with LNPs containing FLuc mRNA (5µg/leg) or saRNA (1µg/leg) and expression was determined using fluorescence and luminescence imaging (IVIS) over 28 days. D-luciferin (150mg/kg) was injected IP 15-20 minutes before luminescence imaging. To determine the immunogenicity of the saRNA-LNPs, female BALB/c mice (n=5) were immunized by IM injection on day 0 with LNPs encapsulating 1µg nCoV saRNA and boosted at day 28. IgG levels in serum on day 27 and day 42 were measured by ELISA.



GenVoy-ILM is an Effective *In Vivo* Delivery Vehicle for mRNA and saRNA

LNP Delivery (A): LNP delivery and clearance from the muscle is comparable for both GenVoy-ILM and MC3 LNPs over 28 days. **Luciferase Expression (B):** Luciferase expression from FLuc mRNA and saRNA over 28 days, delivered using GenVoy-ILM and MC3 LNPs. N=5 ± SD



Conclusion

- GenVoy-ILM is a commercially available and easy-to-use LNP reagent mix that enables the rapid preclinical development of RNA vaccines.
- GenVoy-ILM LNPs can be used to encapsulate both conventional mRNA and self-amplifying RNA while retaining key critical quality attributes for an RNA vaccine.
- GenVoy-ILM LNPs are an effective *in vitro* and *in vivo* delivery vehicle for both mRNA and saRNA.
- GenVoy-ILM LNPs encapsulating SARS-CoV-2 saRNA induce an immune response in mice, which highlights its utility as a vehicle for screening RNA in preclinical development.

Acknowledgments

Funding supported by SMART:SCOTLAND