

A Novel RNA Lipid Nanoparticle Platform: Gene-edited CAR T Cells for Off-the-Shelf Cancer Therapy

R. Geczy^{a,b}, A. Balgi^a, S. Park^a, R. Zhao^a, C. Webb^a, E. Watt^a, M. Fujisawa^a, N. Jain^a, A. Zhang^a, A. Geall^a, A. Thomas^a, and S. Clarke^{a,c}

^a Affiliation: Precision NanoSystems ULC, 55 W Kent Ave N, Vancouver, BC V6P 6T7


^b Presenting author email: reka_geczy@pall.com | ^c corresponding author email: samuel_clarke@pall.com

Contact us at:
Precision NanoSystems ULC, Vancouver, BC, Canada
info@precision-nano.com

Introduction

- Engineering T cells to express chimeric antigen receptor (CAR) redirects these cells to target tumours, making this a promising cell-based cancer therapy [1].
- To engineer universal CAR T cells for allogeneic cell therapy, gene editing can be used to remove risks associated with graft-versus-host disease [2,3].
- However, conventional delivery methods such as viral vectors and electroporation have performance and manufacturing shortcomings that limit their utility for the multi-step engineering of T cells [4].
- Lipid nanoparticles (LNPs) solve these shortcomings, using endogenous uptake pathways to deliver therapeutic RNAs, making this technology gentle to cells. LNPs are fully synthetic and can be rapidly scaled-up using microfluidics [5,6].
- Here we report on the use of a novel LNP reagent in an optimized protocol to achieve successful multi-step engineering of primary T cells with high efficiency and high cell viability.

LNP-based Off-the-Shelf T Cell Therapy




LNP-based delivery of nuclease and therapeutic mRNAs

- Gene editing with:
- CRISPR-Cas9
 - TALENs/Zinc fingers
- Combined with:
- CAR / TCR mRNA

Material and Methods

Multi-step Engineering for Gene-Edited CAR T Cells



Materials

Lipid nanoparticle reagent, GenVoy-ILM™ T Cell Kit for mRNA, Precision NanoSystems ULC, 1000683 and 1001144 • CleanCap® Cas9 mRNA (5moU), Trilink, L-7206 • CleanCap® CD19-CAR mRNA (wt), Trilink, custom order • sgRNAs Alt-R® CRISPR-Cas9 sgRNA, Integrated DNA Technologies, custom targets • Quant-iT™ RiboGreen® RNA Assay Kit, Thermo Fisher Scientific, R11490 • Cryopreserved T cells, STEMCELL Technologies Inc., 70024 • ImmunoCult™ Human CD3/CD28/CD2 T Cell Activator, STEMCELL Technologies Inc., 10970 • ImmunoCult™-XF T Cell Expansion Medium, STEMCELL Technologies Inc., 10981

Equipment

NanoAssemblr Spark (NIS0003) and Ignite (NIS0001), Precision NanoSystems ULC • BioTek™ Synergy™ H1 plate reader • Cytoskeleton V3-B3-R0, Beckman Coulter, C09747

Results

Expression on CD19 CAR with High Efficiency and High Cell Viability




Figure 1. A) Schematic illustration of the custom CAR mRNA construct. B) CAR transfection efficacy 24 h post LNP addition in primary T cells. Assessed using flow cytometry with n=2 sample replicates. C) Surface expression profile of CAR analyzed by flow cytometry. D) Cell viability normalized to untreated (UT) cells. Significance evaluated using one-way ANOVA.

Knock Out of TCR and CD3 in a Dose-dependent Manner




Figure 2. A) Percent TCR (red) or CD3 (orange) knockout of TCR sgRNA treated samples at the indicated 0.2 – 4.8 µg RNA / million cells dose. B) Surface expression of TCR as detected by flow cytometry at the indicated RNA dose. C) Percent cell viability normalized to untreated (UT) for LNP treated samples at the indicated doses. Cas9 mRNA and sgRNA were co-encapsulated at a 1:1 wt ratio, 10 µg total RNAs used per LNP sample.

Knock Out of Both CD52 and TCR in a Multiplex Gene Editing Approach




Figure 3. A) Schematic illustration of the experiment. For multiplexing, two LNP samples were prepared, one containing TCR sgRNA and with one CD52 sgRNA. B) Percent TCR and CD52 double negative cells at a 3.2 µg/million cells RNA dose. C) Corresponding quadrant analysis of CD52 and TCR levels in the untreated (UT) or LNP treated T cell population.

Conclusions

- The presented LNP reagent, GenVoy-ILM T Cell Kit for mRNA, is versatile for gene expression, gene editing, or the combination of both in a multi-step approach.
- LNPs show exceptional performance with $81.5 \pm 4.3\%$ CD19 CAR expression and $74.7 \pm 5.1\%$ TCR gene knockout in primary T cells. LNPs are gentle and maintain $>90\%$ cell viability.
- LNPs are easily scalable using PNI's microfluidic technology, where both the discovery and preclinical LNPs show excellent performance.

Multi-step Gene-edited CAR T Cells are Highly Functional for Cancer Therapy




Figure 4. A) TCR KO on Day 7 when treated with TCR sgRNA and Cas9 mRNA at a dose of 3.2 µg RNA/million cells. Starting sample was purified by TCR negative selection prior to downstream applications. B) Percent CD19-CAR expression 24 h after treatment with CAR mRNA LNPs at 3.2 µg RNA/million cells, in the untreated or TCR KO populations. Lower CAR expression is attributed to a lower activation state after TCR knockout. C) Percent cell viability normalized to UT after CAR expression or CAR expression and TCR knockout. D) Specific lysis of SUP-B15 cells by UT, TCR+/CAR19+, or the gene edited TCR-/CAR19+ T cells at the indicated effector (T cell, E) to target (B cell, T) ratios, E:T. Significance evaluated using t-tests (Fig. B) or one-way ANOVA with multiple comparisons (Fig. D).

Scalable Microfluidic Technology Allows for Discovery to Pre-clinical Translation




Figure 5. A) Seamless scale-up allows for the production of discovery scale LNPs at 10 µg RNA, to pre-clinical scale at 4 mg RNA. B) Nanoparticle Tracking Analysis (NTA) shows agreeable size distribution between the two scales. C) A CD19-CAR transfection efficiency 24 h post mRNA LNP addition for both scales. D) Single target (T cell receptor, TCR) knockout, and E) double target (TCR and CD52) knockout through sgRNA and Cas9 mRNA delivery. F) Functional killing of CD19+ B cells (SUP-B15) in 16 h co-culture experiment for both scales. For all: a dose of 3.2 µg RNA/million cells was applied to human primary T cells. Average gene expression and gene knockout shows at least n=8 LNPs and n=2 donors. Functional killing shows n=2 LNPs and n=2 donors. Error bars represent standard deviation with significance evaluated using t-tests among selected groups.

References

- [1] Feins, S., et al., An introduction to chimeric antigen receptor (CAR) T-cell immunotherapy for human cancer. *Am J Hematol*, 2019. 94(S1): p. S3-S9. [2] Rafiq, S., C.S. Hackett, and R.J. Brentjens, Engineering strategies to overcome the current roadblocks in CAR T cell therapy. *Nat Rev Clin Oncol*, 2020. 17(3): p. 147-167. [3] Zhao, J., et al., Universal CARs, universal T cells, and universal CAR T cells. *J Hematol Oncol*, 2018. 11(1): p. 132. [4] Atsavapranee, E.S., M.M. Billingsley, and M.J. Mitchell, Delivery technologies for T cell gene editing: Applications in cancer immunotherapy. *EBioMedicine*, 2021. 67: p. 103354. [5] Absalon, J., K. Koury, and W.C. Gruber, Safety and Efficacy of the BNT162b2 mRNA Covid-19 Vaccine. *Reply. N Engl J Med*, 2021. 384(16): p. 1578. [6] Baden, L.R., et al., Efficacy and Safety of the mRNA-1273 SARS-CoV-2 Vaccine. *N Engl J Med*, 2021. 384(5): p. 403-416.

Acknowledgments

We thank the contributions of Rebecca De Souza, Andrew Brown and Helena Son during the early development of the T Cell Kit. We further thank the guidance of Vicki Stronge throughout the development process, and Suraj Abraham for assistance in the design of the CAR construct. Additionally, Ghinwa Darwish for conducting the NTA measurements.