

# Well Characterized Lipid Nanoparticle Library Accelerates Development of Next Generation Genomic Medicine

Serra Gurcan<sup>1</sup> (Presenting Author), Nikita Jain<sup>1</sup> (Corresponding Author), Sedigheh Nazaripour<sup>1</sup>, Zhengyu Chen<sup>1</sup>, Suraj Abraham<sup>1</sup>, Leanna Yee<sup>1</sup>, Sams Sadat<sup>1</sup>, Gayatri Mehar Namala<sup>1</sup>, Sijo Chemmannur<sup>1</sup>, Kobe Tickner<sup>1</sup>, Malathi Anantha<sup>1</sup>, Ruchi Sharma<sup>1</sup>, Srinivas Abbina<sup>1</sup>, Seetalakshmi Thambatti<sup>1</sup>, Vinay Mayya<sup>1</sup>, Emily Soon<sup>1</sup>, Jay Paquette<sup>1</sup>, Anitha Thomas<sup>1</sup>

<sup>1</sup>Precision Nanosystems ULC, Vancouver, Canada


Contact us at:  
Precision NanoSystems, Vancouver, BC, Canada  
info@precision-nano.com

## Introduction

- FDA approval of ONPATRO<sup>®</sup> by Alnylam, Comirnaty<sup>®</sup> by BioNTech/Pfizer and Spikevax<sup>®</sup> by Moderna and the various clinical trials with mRNA-based drugs or vaccines have provided momentum to further develop lipid nanoparticle (LNP) based genetic medicines.
- Ionizable amino lipids are a major constituent of LNPs for delivering nucleic acid therapeutics, and thus ionizable lipids with high encapsulation efficiency, high endosomal release that are non-toxic are essential for efficient clinical translation.
- The scarcity of ionizable lipids that are suitable for development of vaccines, cell and gene therapies continues to be a problem in advancing many potential therapeutic/vaccine candidates to the clinic.



## Methods and Results

### 1. Non-Viral Delivery Systems Enable Cell Therapy



- MC3 – Ionizable lipid, an excipient in FDA approved Onpatro<sup>™</sup>, was used as benchmark lipid for comparative evaluation of potency of proprietary PNI lipids
- Higher GFP MFI and comparable transfection efficiency relative to MC3 was observed with PNI lipids
- Cell viability was >90% as compared to untreated cells




### 2. Proprietary LNPs Enable Erythropoietin Production *in vivo*



- EPO-encoded mRNA-LNP showed significant EPO expression levels at 6h in C57BL/6 mice following i.v administration of 0.5 mg/kg dose
- Post 7 days injection, EPP encoded mRNA LNP treated female C57BL/6 mice demonstrated ~20–40% increase in Hematocrit levels

### 3. Proprietary LNPs Towards Developing Self-Amplifying mRNA (saRNA) Vaccines


How do saRNA-LNP vaccines work?



- The analysis for SARS CoV-2 spike protein specific antibodies, post two weeks boost (day 42), confirmed that many of the PNI lipids showed similar expression compared to clinically demonstrated SM-102 and ALC-0315.


## Objectives

- Demonstrate cell therapy applications for proprietary lipids using GFP encoded mRNA and compare their transfection with clinically approved lipid Dlin-MC3-DMA (MC3) in human primary T cells.
- Demonstrate protein replacement applications for proprietary lipids using EPO encoded mRNA and their potency comparison with clinically approved Dlin-MC3-DMA in mice.
- Display PNI proprietary lipids for vaccine applications using self-amplifying RNA encoding for SARS-CoV-2 spike protein in comparison to SM-102 and ALC-0315 in mice.
- Showcase the PNI proprietary lipids for gene editing applications *in vivo*.
- Illustrate the safety and tolerability of PNI proprietary LNPs in mice.




### 4. Novel Ionizable Lipids for Gene Editing Applications

#### a. TTR CRISPR Gene Editing – Preliminary Data with PNI Delivery System



#### b. Base Editing - PCSK9 Knock Down - Preliminary Data with PNI Delivery system




### 5. Tolerability of Proprietary LNP Administered IV in Mice

72h study CD-1 male mice 6wk old n=4/LNP	0-72 h post dose		72 h post dose			
	Clinical Observations		Necropsy		Liver Histology	
	Moderate Dose	High Dose	Moderate Dose	High Dose	Moderate Dose	High Dose
Control PBS	Normal		Normal		Normal	
SM-102	Mild signs recover w/in 6h	Moderate signs lasted full 72h	Normal	3/4 Normal 1/4 Euthanized @ 1h	Normal	Slight apoptosis & inflammation
ALC-0315	Mild signs recover w/in 72h	Mild to moderate for 6-72h	Normal	3/4 Normal 1/4 Enlarged Spleen	Mild apoptosis & inflammation	Mild apoptosis & inflammation
PNI iL2	Normal	Mild signs recover w/in 6h	Normal	Normal	Slight apoptosis & inflammation	Slight apoptosis & inflammation

**Clinical Observations** were relatively higher, lasted longer, and occurred with greater frequency in the LNPs with SM-102 and ALC-0315, than compared to a matched mg/kg dose of Precision NanoSystems LNP.

**Necropsy** were normal for most groups.

**Liver Histology** showed relatively higher apoptosis and inflammation in the LNP containing ALC-0315 compared to all other groups at matched doses.



## Conclusion

- Non-viral lipid nanoparticle delivery systems show significant promise in the field of genomic medicine.
- Precision NanoSystems has developed a proprietary ionizable lipid library comprising more than 100 lipids with diverse pKa for different applications including cell therapy, protein replacement, gene therapy and RNA vaccine.
- The Precision NanoSystems lipid technology enables the targeted delivery of nucleic acids to specific cells and tissues and can help to accelerate the development of genomic medicines for a wide-range of diseases.



Download