Accelerated development of self-amplifying RNA (saRNA) vaccines using scalable NxGen[™] Microfluidics

Edward Grahame (Presenting Author), Ariel Zhang¹, Sams Sadat¹, Pierrot Harvie¹, Nikita Jain¹, Sijo Chemmannur¹, Leanna Yee¹, Srinivas Abbina¹, Helena Son¹, Cece Nazaripour¹, Suraj Abraham¹, Sita Thampatty¹, Darius Menezes², Alex Shapiro¹, Derrick Stebbing¹, Emily Soon¹, Karen Markwick¹, Aleksei Angell¹, Jay Paquette¹, Anna Blakney³, Robin Shattock³, Andrew Geall¹, Lloyd Jeffs¹, Anitha Thomas¹

¹Precision NanoSystems, Vancouver, BC, Canada ²Cytiva, Vancouver, BC, Canada ³Imperial College London, London, United Kingdom

Introduction

- The COVID-19 pandemic has unveiled the potential of messenger RNA (mRNA) based vaccines as an ideal platform for pandemic response.
- RNA vaccines mimic antigen structure and expression similar to natural infection without causing the disease, but allowing body to produce antibody response against future infection.
- saRNAs have the potential for antigen sparring since lower doses elicit effective immune response compared to non-replicating mRNA vaccines.
- o saRNAs are inherently more fragile than mRNA and are prone to degradation. Typical high shear manufacturing methods affect the potency.
- Herein, we showcased the utility of microfluidics to enable low shear, rapid screening of preclinical candidates and the swift advancement to GMP-enabling studies.
- The LNP downstream processing was also optimized using the Pall delta TFF cassette. These Pall cassette showed high flux as compared to competitor product.

How do saRNA-LNP vaccines work?



Printer

Objectives

- Production of Self-Amplifying mRNA (saRNA) SARS-CoV-2 Vaccine.
- Evaluate cellular and humoral responses of the vaccine candidates.
- Evaluate the efficacy of the vaccine candidates in a SARS-CoV-2 Hamster challenge model.

Methods

- SARS-COV-2 full length spike protein encoded self-amplifying RNAs were encapsulated in Lipid nanoparticles using low shear NxGen[™] microfluidics platform.
- Two RNA-vaccine candidate formulations (LNP-1 and LNP-2) manufactured using the Precision NanoSystems Manufacturing Platform (Ignite[™], Blaze[™], NanoAssemblr[®] GMP) were assessed for scalability.



Booth #40

Contact us at: Precision NanoSystems, Vancouver, BC, Canada info@precision-nano.com

saRNA is enzymatically synthesized in vitro from a DNA template

saRNA is delivered in a LNP formulation, to cytoplasm of host cells







•_• • •

Results

1. Vaccine candidates manufactured using NxGen[™] Microfluidics show effective Humoral and T Cell Responses





CCL4 (MIP-1β) CCL5(RANTES) IFNg



2. Two LNP compositions, LNP-1 and LNP-2 were selected for scale up based on *in vitro* and *in vivo* activity



Manufacturing Process

	Ignite™	Blaze™	GMP					
NxGen™ Mixers	NxGen, NxGen w/in-line dilution	NxGen 400, NxGen 500	NxGen 500					
Organic Phase	Lipid in Ethanol							
Aqueous Phase	RNA in aqueous							
Total Micromixing Maximum Volume	up to 20mL undiluted	up to 10L undiluted	up to 50L undiluted					
Flow Rate Ratio [Org : Aq]	3:1							
Total Flow Rate	12 mL/min	115 mL/min	200 mL/min					
In-line Dilution Ratio (Buffer : Micromix volume)	3:1	3:1	3:1					
Downstream Processing	Dialysis Cassette	TFF	TFF					

8. Vaccine candidates protected hamsters in a SARS-CoV-2 challenge study





Tissue viral RNA

Cranial

Vaccinated hamsters had lower SARS-CoV-2 viral loads in Nasal wash and Respiratory tissues



10³ Study day A = Vehicle control B = Vaccine Candidate 1 C = Vaccine Candidate 2

Vaccinated hamsters had higher levels o

SARS-CoV-2 specific IgG antibody and

neutralizing antibody titers

9. Immune response of vaccine candidates in NHP

3. Self-amplifying RNA-LNP have equivalent Size, PDI and encapsulation across scales (Ignite-Blaze-GMP)



SARS-CoV-2 self-amplifying RNA-LNP made with Precision NanoSystems proprietary ionizable lipid had similar size (\sim 60 nm), polydispersity (<0.2) and encapsulation efficiency (>90%) across all scales tested with two different LNPs

4. Scalability Across Platforms was achieved for LNP-1 and LNP-2

Simplified Scale-Up of mRNA-LNP Using NxGen[™]





5. LNP composition was selected based on stability and robustness during TFF process

• Both LNP-1 and LNP-2 efficiently induced SARS-CoV-2 specific IgG response in mice

• As observed *in vitro* LNP-1 showed slightly higher IgG response as compared to LNP-2

• Both LNP-1 and LNP-2 generated neutralizing antibodies against the SARS-CoV-2 virus



• At the 50 mg scale LNP-1 showed particle size increase during the downstream processing (TFF)

• LNP-2 remained size stable during the downstream processing step and was selected as the lead for GLP tox and GMP manufacturing



- Vaccine candidates showed high SARS-CoV-2 specific IgG titres in NHPs
- Significant difference in IgG expression at Day 42 (post prime) was observed compared to pre-bleeds and the response was amplified post booster (Day 70)

10. Long term storage stability of vaccine candidate and morphology



- Vaccine candidate was stable up to 12-months for similar immune response in animals
- All the LNPs throughout this study were found in the range of 70-75 nm with encapsulation efficiency more than 95%

TEM image of Precision NanoSystems LN vaccine candidate at higher magnificatior 73.000x) with scale bar of 200 nm.

11. Precision NanoSystems' Lyophilized saRNA-LNP Vaccine Candidate Retains Activity

• Similar particle characteristics and in vitro potency following lyophilization cycle

	Lyophilization	Lyophilization
Size (d.nm)	71	89
PDI	0.074	0.09
% saRNA encapsulation	97	95
EC50 (ug/mL)	0.063	0.057



BHK 570 cells transfected in a 96-well plate with SARS-CoV-2 saRNA LNP in a dose response manner from 1 to 0.00049 ug/mL

Pre Lyophilization

Post-Lyophilizatior

Conclusion

• LNP based Vaccine candidates manufactured using Precision NanoSystems NanoAssemblr® platform showcased effective cellular and humoral immune response.

6. LNP-2 Downstream processing optimization by flux excursion



Precision NanoSystems COVID-19 vaccine candidate in 3% ethanol was recirculated into the TFF product for 5 min at various feed flow rate prior to TMP increase in order to generate these flux excursion curve

%saRNA Encapsulation PDI **TFF Product** T=0 Final T=0 Final T=0 Final Pall T-series 69 94 0.16 0.19 99 99 Centramate™cassette, Delta part# DC030T01 Cassette X 300 KD 71 0.19 0.25 99 99 62 MWCO Cassette Y 100 KD 105 0.17 0.21 99 65 98 MWCO Hollow fiber X 300 KD 69 74 0.11 0.14 99 99 MWCO

Flux excursion demonstrated that the Pall T-series Centramate[™]cassette, Delta (DC030T01) was more effective in terms of permeate flux for LNP diafiltration as compared to competitor TFF product

The Pall cassette can operate at higher TMP thus leading to a superior permeate flux without significant changes in LNP characteristics

7. The Pall T-series Centramate[™] cassette, Delta performs better than a hollow fiber product at the 50 mg scale

LNP Processing Time	PALL Delta cassette 1000cm ² MWCO 30KD	Hollow Fiber 1600cm ² MWCO 500KD		Size (nm)		PDI		% Encap.	
				Cassette	Hollow Fiber	Cassette	Hollow Fiber	Cassette	Hollow Fiber
TFF conc step (6.0X – 6.5X) ~33mir	~33min	~20 min	8X dilution	65	59	0.172	0.104	99	99
			TFF concentrated	73	56	0.256	0.091	99	98
Diafiltration in storage buffer (4DV)	~37 min	6.0X – 6.5X	/5	50	0.230	0.051	55	50	
		TFF dilution buffer	83	69	0.278	0.098	98	97	
Final conc step ~5min	~6 min								
			TFF final collected	92	68	0.343	0.083	98	97
Total time	1h	1h 1h 07min	Sterile filtered	77	69	0.17	0.076	98	97

• The total TFF time on a Pall delta 30 KD is 7 min faster with a smaller filter as compared to processing time on a larger hollow filter

- Both TFF filters lead to comparable particles characteristics
- The LNP from the Pall cassette were sterile filtered on a Pall ECV and the LNP from the hollow fiber TFF were sterile filtered on a Sartopore 2

- SARS-CoV-2 self-amplifying RNA-LNP made with Precision NanoSystems proprietary ionizable lipid had similar Critical Quality Attributes (CQAs) such as size (~70 nm), polydispersity (<0.2) and encapsulation efficiency (>90%) across all scales tested with two different LNP compositions.
- Precision NanoSystems COVID-19 vaccine candidate protected hamsters in a SARS-CoV-2 challenge study and showed high SARS-CoV-2 specific IgG titres in NHPs.
- SARS-CoV-2 self-amplifying RNA LNP Vaccine candidates showed long term stability and is compatible with a lyophilization process.
- Downstream processing time and particle stability during large scale TFF should be considered as a critical parameter during scale up.
- Pall product showed faster downstream processing time as compared to other membrane tested.

Acknowledgments

We kindly acknowledge the work of various team members in Robin Shattock's Lab, Imperial College London, and Delivery, RNA Services, QC, and Clinical manufacturing teams in Precision NanoSystems. We also thank the generous funding support from Strategic Initiative Fund, Canada, without which this work would not have happened.





